

Product Datasheet

Human ILK-1 ELISA Kit (orb548818)

Catalog Number	orb548818
Category	Assays and Kits
Description	Human Integrin Linked Kinase 1 (ILK-1) Sandwich ELISA Kit
Reactivity	Human
Form/Appearance	Traditional
Assay Type	Sandwich
Assay Time	4.5 hours
Range	78.1-5000pg/mL
Sensitivity	33.5pg/mL
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)