

Product Datasheet

Human IL22Ra2 ELISA Kit (Ready to Use) (orb548831)

Catalog Number	orb548831
Category	Assays and Kits
Description	Human Interleukin 22 Receptor Alpha 2 (IL22Ra2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient and simple assays.
Reactivity	Human
Form/Appearance	Ready to Use
UniProt ID	Q969J5
Assay Type	Sandwich
Assay Time	3 hours
Range	0.312-20ng/mL
Sensitivity	0.127ng/mL
Sample Types	serum, plasma, tissue homogenates and other biological fluids
Storage	Store at -20°C
Note	For research use only
Entrez	116379

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)