

Product Datasheet

Other VB6 ELISA Kit (Ready to Use) (orb549600)

| | |
|--------------------------|--|
| Catalog Number | orb549600 |
| Description | General Vitamin B6 (VB6) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient and simple assays. |
| Reactivity | All |
| Range | 3.12-200ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Sample Type: serum, plasma, tissue homogenates and other biological fluids; competitive inhibition: Experimental Method |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Assay Time | 3 hours |
| Sensitivity | 1.17ng/mL |
| Expiration Date | Please enquire. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)