

Product Datasheet

Human AP36 ELISA Kit (Ready to Use) (orb549858)

Catalog Number	orb549858
Category	Assays and Kits
Description	Human Apelin 36 (AP36) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient and simple assays.
Reactivity	Human
Form/Appearance	Ready to Use
Assay Type	Competitive
Assay Time	3 hours
Range	46.8-3000 pg/mL
Sensitivity	12.58 pg/mL
Sample Types	serum, plasma and other biological fluids
Storage	Store at -20°C
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)