

## Product Datasheet

### Human ABCG1 ELISA Kit (Ready to Use) (orb549939)

|                        |                                                                                                                                                     |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>  | orb549939                                                                                                                                           |
| <b>Category</b>        | Assays and Kits                                                                                                                                     |
| <b>Description</b>     | Human ATP Binding Cassette Transporter G1 (ABCG1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient and simple assays. |
| <b>Reactivity</b>      | Human                                                                                                                                               |
| <b>Form/Appearance</b> | Ready to Use                                                                                                                                        |
| <b>UniProt ID</b>      | <b>P45844</b>                                                                                                                                       |
| <b>Assay Type</b>      | Sandwich                                                                                                                                            |
| <b>Assay Time</b>      | 3 hours                                                                                                                                             |
| <b>Range</b>           | 0.156-10ng/mL                                                                                                                                       |
| <b>Sensitivity</b>     | 0.061ng/mL                                                                                                                                          |
| <b>Sample Types</b>    | tissue homogenates, cell lysates and other biological fluids                                                                                        |
| <b>Storage</b>         | Store at -20°C                                                                                                                                      |
| <b>Note</b>            | For research use only                                                                                                                               |
| <b>Entrez</b>          | <b>9619</b>                                                                                                                                         |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)