

Product Datasheet

Human LDHA ELISA Kit (Ready to Use) (orb551325)

Catalog Number	orb551325
Category	Assays and Kits
Description	Human Lactate Dehydrogenase A (LDHA) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient and simple assays.
Reactivity	Human
Form/Appearance	Ready to Use
UniProt ID	P00338
Assay Type	Sandwich
Assay Time	3 hours
Range	0.781-50ng/mL
Sensitivity	0.32ng/mL
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Entrez	3939

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)