

## Product Datasheet

### Human MHCA ELISA Kit (Ready to Use) (orb551446)

|                        |                                                                                                                                                          |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>  | orb551446                                                                                                                                                |
| <b>Category</b>        | Assays and Kits                                                                                                                                          |
| <b>Description</b>     | Human Major Histocompatibility Complex Class I A (MHCA) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient and simple assays. |
| <b>Reactivity</b>      | Human                                                                                                                                                    |
| <b>Form/Appearance</b> | Ready to Use                                                                                                                                             |
| <b>UniProt ID</b>      | <b>P30443</b>                                                                                                                                            |
| <b>Assay Type</b>      | Sandwich                                                                                                                                                 |
| <b>Assay Time</b>      | 3 hours                                                                                                                                                  |
| <b>Range</b>           | 0.312-20ng/mL                                                                                                                                            |
| <b>Sensitivity</b>     | 0.126ng/mL                                                                                                                                               |
| <b>Sample Types</b>    | tissue homogenates and other biological fluids.                                                                                                          |
| <b>Storage</b>         | All reagents should be stored at -20°C upon receipt                                                                                                      |
| <b>Note</b>            | For research use only                                                                                                                                    |
| <b>Entrez</b>          | <b>3105</b>                                                                                                                                              |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)