

Product Datasheet

Human SUMF1 ELISA Kit (Ready to Use) (orb552486)

Catalog Number	orb552486
Category	Assays and Kits
Description	Human Sulfatase Modifying Factor 1 (SUMF1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient and simple assays.
Reactivity	Human
Form/Appearance	Ready to Use
UniProt ID	Q8NBK3
Assay Type	Sandwich
Assay Time	3 hours
Range	0.156-10ng/mL
Sensitivity	0.064ng/mL
Sample Types	tissue homogenates and other biological fluids.
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Entrez	285362
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)