

Product Datasheet

Mouse CPT1A ELISA Kit (Ready to Use) (orb553097)

Catalog Number	orb553097
Category	Assays and Kits
Description	Mouse Carnitine Palmitoyltransferase 1A, Liver (CPT1A) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient and simple assays.
Reactivity	Mouse
Form/Appearance	Ready to Use
UniProt ID	P97742
Assay Type	Sandwich
Assay Time	3 hours
Range	0.156-10ng/mL
Sensitivity	0.055ng/mL
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Entrez	12894

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)