

Product Datasheet

Mouse IL1F9 ELISA Kit (Ready to Use) (orb553581)

Catalog Number	orb553581
Category	Assays and Kits
Description	Mouse Interleukin 1 Family, Member 9 (IL1F9) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient and simple assays.
Reactivity	Mouse
Form/Appearance	Ready to Use
UniProt ID	Q3U0P4
Assay Type	Sandwich
Assay Time	3 hours
Range	7.81-500pg/mL
Sensitivity	2.6pg/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Entrez	215257

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)