

Product Datasheet

Mouse PGC ELISA Kit (Ready to Use) (orb553883)

Catalog Number	orb553883
Category	Assays and Kits
Description	Mouse Pepsinogen C (PGC) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient and simple assays.
Reactivity	Mouse
Form/Appearance	Ready to Use
UniProt ID	Q9D7R7
Assay Type	Sandwich
Assay Time	3 hours
Range	0.156-10ng/mL
Sensitivity	0.063ng/mL
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Entrez	109820

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)