

Product Datasheet

Rat GMCSF ELISA Kit (Ready to Use) (orb554813)

Catalog Number	orb554813
Category	Assays and Kits
Description	Rat Colony Stimulating Factor 2, Granulocyte Macrophage (GMCSF) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient and simple assays.
Reactivity	Rat
Form/Appearance	Ready to Use
UniProt ID	P48750
Assay Type	Sandwich
Assay Time	3 hours
Range	15.6-1000pg/mL
Sensitivity	6.5pg/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Entrez	116630

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)