

Product Datasheet

Rat ITaC ELISA Kit (Ready to Use) (orb555123)

Catalog Number	orb555123
Category	Assays and Kits
Description	Rat Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient and simple assays.
Reactivity	Rat
Form/Appearance	Ready to Use
UniProt ID	Q7TNL0
Assay Type	Sandwich
Assay Time	3 hours
Range	46.8-3000pg/mL
Sensitivity	16.13pg/mL
Sample Types	serum, plasma and other biological fluids
Storage	Store at -20°C
Note	For research use only
Entrez	305236

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)