

Product Datasheet

Rabbit IgG / PE (orb555865)

Catalog Number orb555865

Category Proteins

Description This product is Rabbit IgG / PE, also known as Rabbit Immunoglobulin G/PE

Isotype Control. It is provided in liquid form.


Form/Appearance Liquid

Buffer/Preservatives 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% glycerol.

Storage Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note For research use only

Expiration Date 12 months from date of receipt.



Blank control: SH-SY5Y. Isotype control: Rabbit IgG-PE (orb555865), dilution: 1 ug/Test, Primary Antibody (green line): Rabbit Anti-ChAT antibody, dilution: 1 ug/Test, Protocol The cells were fixed with 4% PFA (10 min at room temperature) and then permeabilized with 90% ice-cold methanol for 20 min at -20°C. The cells were then incubated in 5% BSA to block non-specific protein-protein interactions for 30 min at room temperature. Cells stained with Primary Antibody for 30 min at room temperature. Isotype control antibody (orange line) was rabbit IgG. Acquisition of 20000 events was performed.