


Product Datasheet

Mouse 8-iso-PGF2alpha ELISA Kit (orb566240)

Catalog Number	orb566240
Category	Assays and Kits
Description	Mouse 8-iso-PGF2alpha(8-isoprostane) ELISA Kit
Target	8-iso-PGF2 α
Reactivity	Mouse
Assay Type	Competitive
Assay Time	2 hours
Range	15.625-1000pg/ml
Sensitivity	9.375pg/ml
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Storage	2 - 8 °C for 12 months
Note	For research use only



Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)