


Product Datasheet

Rat Hspa8 ELISA Kit (orb566940)

Catalog Number	orb566940
Description	Rat Hspa8 ELISA Kit
Reactivity	Rat
Range	0.625-40ng/ml
Target	Hspa8
Storage	2 - 8 °C for 12 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Assay Time	4 hours
Uniprot ID	P63018
Sensitivity	0.375ng/ml
Expiration Date	Please enquire.



Rat Hspa8(Heat shock cognate 71 kDa protein) ELISA Kit

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)