


Product Datasheet

Rat 5-hydroxytryptamine ELISA Kit (orb567754)

Description	Rat 5-hydroxytryptamine ELISA Kit
Reactivity	Rat
Range	1.563-100ng/ml
Target	5-hydroxytryptamine
Storage	2 - 8 °C for 12 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Assay Time	2 hours
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Rat 5-hydroxytryptamine ELISA Kit

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)