

Product Datasheet

Saphir LAMP GreenMaster (orb572666)

Catalog Number	orb572666
Category	Tools
Description	Saphir Bst GreenMaster is a complete 2x conc. master mix for isothermal amplification of DNA. The mix is based on a genetically optimized Bst polymerase that allows rapid and specific amplification of DNA at constant temperature (60 to 65 °C). The enzyme shows high strand displacement activity and generates an amplification factor of up to 10 ⁹ which is comparable to approx. 30 cycles in a PCR assay. LAMP technique allows detection of a target gene within 10-30 minutes.
Form/Appearance	liquid
Concentration	2x conc
Application notes	Spectroscopic Propertie: λ_{exc} 495 nm (SYBRGreen bound to DNA), λ_{em} 520 nm (SYBRGreen bound to DNA).
Storage	store at -20 °C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)