

## Product Datasheet

### RARB Rabbit Polyclonal Antibody (orb573746)

|                     |                                                                                                                                         |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Catalog Number      | orb573746                                                                                                                               |
| Description         | Rabbit polyclonal antibody to RARB                                                                                                      |
| Species/Host        | Rabbit                                                                                                                                  |
| Reactivity          | Human                                                                                                                                   |
| Conjugation         | Unconjugated                                                                                                                            |
| Tested Applications | WB                                                                                                                                      |
| Immunogen           | The immunogen is a synthetic peptide directed towards the N-terminal region of Human RARB                                               |
| Target              | RARB                                                                                                                                    |
| Preservatives       | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.                                       |
| Form/Appearance     | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.                                       |
| Concentration       | 0.5 mg/ml                                                                                                                               |
| Storage             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note                | For research use only                                                                                                                   |
| Protein Sequence    | Synthetic peptide located within the following region:<br>EKALKACFSGLTQTEWQHRHTAQSIETQSTSSEELVPSPPSPLPPPRVYK                            |
| Clonality           | Polyclonal                                                                                                                              |
| MW                  | 49kDa                                                                                                                                   |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.


68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Uniprot ID**

**P10826**

**Expiration Date**

12 months from date of receipt.



Host: Rabbit  
Target Name: RARB  
Sample Tissue: Uterus Tumor Lysate  
Antibody Dilution: 1.0µg/ml

Sample Type: Uterus Tumor lysates, Antibody dilution: 1.0 ug/ml.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)