

## Product Datasheet

### PAWR Rabbit Polyclonal Antibody (orb574032)


|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb574032                                                                                                                               |
| <b>Description</b>         | Rabbit polyclonal antibody to WT1                                                                                                       |
| <b>Species/Host</b>        | Rabbit                                                                                                                                  |
| <b>Reactivity</b>          | Human, Rat                                                                                                                              |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | IHC, WB                                                                                                                                 |
| <b>Immunogen</b>           | The immunogen is a synthetic peptide directed towards the middle region of human PAWR                                                   |
| <b>Target</b>              | PAWR                                                                                                                                    |
| <b>Preservatives</b>       | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.                                       |
| <b>Form/Appearance</b>     | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.                                       |
| <b>Concentration</b>       | 0.5 mg/ml                                                                                                                               |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Protein Sequence</b>    | Synthetic peptide located within the following region:<br>VNIPAAECLDEYEDDEAGQKERKREDAITQQNTIQNEAVNLLDPGSSYLL                            |
| <b>Clonality</b>           | Polyclonal                                                                                                                              |
| <b>MW</b>                  | 37kDa                                                                                                                                   |

#### Biorbyt Ltd.


7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.


68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Uniprot ID** **O75796****NCBI** **NP\_002574****Expiration Date** 12 months from date of receipt.

Human Muscle



WB Suggested Anti-PAWR Antibody, Positive Control: Lane 1: 30 ug rat striatum homogenate, Primary Antibody Dilution: 1:1000, Secondary Antibody: Anti rabbit-HRP, Secondary Antibody Dilution: 1:10000.



WB Suggested Anti-PAWR Antibody Titration: 0.2-1 ug/ml, Positive Control: Jurkat cell lysate.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)