

Product Datasheet

ELP4 Rabbit Polyclonal Antibody (orb576278)

| | |
|-----------------------------|------------------------------------------------------------------------------------------------------------|
| Catalog Number | orb576278 |
| Category | Antibodies |
| Description | Rabbit polyclonal antibody to ELP4 |
| Target | ELP4 |
| Clonality | Polyclonal |
| Species/Host | Rabbit |
| Conjugation | Unconjugated |
| Reactivity | Mouse |
| Predicted Reactivity | Bovine, Canine, Equine, Human, Porcine, Rabbit, Rat |
| Form/Appearance | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. |
| Concentration | 1.0 mg/ml |
| Buffer/Preservatives | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. |
| Purification | Protein A purified |
| Immunogen | The immunogen is a synthetic peptide directed towards the middle region of mouse ELP4 |
| Protein Sequence | Synthetic peptide located within the following region: QKSVYAEGFDGANPQKKQKNILRIGIQNLGSPWGDICCKENCDNNHRL |
| UniProt ID | Q9ER73 |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| MW | 46kDa |
| Tested applications | WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| NCBI | NP_076365 |
| Expiration Date | 12 months from date of receipt. |



WB Suggested Anti-ELP4 Antibody Titration: 2.5 ug/ml, ELISA
Titer: 1:1562500, Positive Control: NIH/3T3 cell lysate.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)