

Product Datasheet

OSMR Rabbit Polyclonal Antibody (orb579776)

Catalog Number	orb579776
Category	Antibodies
Description	Rabbit polyclonal antibody to OSMR
Target	OSMR
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human
Predicted Reactivity	Canine, Equine, Guinea pig, Rabbit
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Buffer/Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification	Affinity Purified
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human OSMR
Protein Sequence	Synthetic peptide located within the following region: LLEKKTGYSQELAPSDNPHVLVDTLTSHSFTLSWKDYSTESQPGFIQGYH
UniProt ID	Q99650

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

MW 110kDa


Tested applications WB

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


Note For research use only

NCBI [NP_003990](#)

Expiration Date 12 months from date of receipt.



Positive control (+): Human Lung Tumor (T-LU), Negative control (-): Human Liver Tumor (T-LI), Antibody concentration: 3 ug/ml.



WB Suggested Anti-OSMR Antibody Titration: 0.2-1 ug/ml, ELISA Titer: 1:312500, Positive Control: 721_B cell lysate.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558