

Product Datasheet

NPTX2 Rabbit Polyclonal Antibody (orb583129)

Catalog Number	orb583129
Category	Antibodies
Description	Rabbit polyclonal antibody to NPTX2
Target	NPTX2
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human
Predicted Reactivity	Bovine, Canine, Equine, Guinea pig, Mouse, Porcine, Rabbit, Rat, Zebrafish
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Buffer/Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification	Affinity Purified
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human NPTX2
Protein Sequence	Synthetic peptide located within the following region: ELEDEKSLLHNESAHRQKTESTLNALLQRVTELERGNSAFKSPDAFKVS
UniProt ID	P47972

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

MW 47kDa**Tested applications** WB**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.**Note** For research use only**NCBI** **NP_002514**WB Suggested Anti-NPTX2 Antibody Titration: 0.2-1 ug/ml,
ELISA Titer: 1:312500, Positive Control: Transfected 293T.**Biorbyt Ltd.**7 Signet Court, Swann Road
Cambridge
CB5 8LA
United KingdomEmail: info@biorbyt.com, support@biorbyt.comPhone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)**Biorbyt LLC**68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United StatesEmail: info@biorbyt.com, support@biorbyt.comPhone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)