

Product Datasheet

PTX3 Rabbit Polyclonal Antibody (orb583189)

Description	Rabbit polyclonal antibody to PTX3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human PTX3
Target	PTX3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: RMLLQATDDVLRGELQRLREELGRLESLARPCAPGAPAEARLTSALDEL
Clonality	Polyclonal
MW	42kDa
Uniprot ID	P26022

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)


NCBI**NP_002843****Expiration Date**

12 months from date of receipt.




Host: Rabbit
Target Name: PTX3
Sample Type: HeLa Cell Lysate
Antibody Dilution: 1.0µg/ml

Sample Tissue: Human HeLa Whole Cell, Antibody dilution: 1 ug/ml.



Host: Rabbit
Target Name: PTX3
Sample Type: Jurkat Cell Lysate
Antibody Dilution: 1.0µg/ml

Sample Tissue: Human Jurkat Whole Cell, Antibody dilution: 1 ug/ml.

PTX3

Host: Rabbit
Target Name: PTX3
Sample Type: MCF7
Antibody Dilution: 1ug/ml


Sample Tissue: Human MCF7, Antibody dilution: 1.0 ug/ml.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Sample Type: HepG2 Whole Cell, Lane A: Primary Antibody, Lane B: Primary Antibody + Blocking Peptide, Primary Antibody Concentration: 1 µg/ml, Peptide Concentration: 5 µg/ml, Lysate Quantity: 25 µg/lane, Gel Concentration: 0.12%.



prostate



WB Suggested Anti-PTX3 Antibody Titration: 0.2-1 µg/ml, ELISA Titer: 1:312500, Positive Control: MCF7 cell lysate.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)