

## Product Datasheet

### MCM6 Antibody (orb607526)

|                            |                                                                                                                                                                                  |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb607526                                                                                                                                                                        |
| <b>Category</b>            | Antibodies                                                                                                                                                                       |
| <b>Description</b>         | Mouse monoclonal antibody to MCM6 (Proliferation Marker)                                                                                                                         |
| <b>Target</b>              | MCM6                                                                                                                                                                             |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                                       |
| <b>Species/Host</b>        | Mouse                                                                                                                                                                            |
| <b>Isotype</b>             | IgG2b                                                                                                                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                     |
| <b>Reactivity</b>          | Human                                                                                                                                                                            |
| <b>Form/Appearance</b>     | 200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml. |
| <b>Purification</b>        | Protein A                                                                                                                                                                        |
| <b>Immunogen</b>           | Recombinant fragment (around aa 228-368) of human MCM6 protein (exact sequence is proprietary)                                                                                   |
| <b>UniProt ID</b>          | <b>Q14566</b>                                                                                                                                                                    |
| <b>MW</b>                  | 105kDa                                                                                                                                                                           |
| <b>Tested applications</b> | IHC, WB                                                                                                                                                                          |
| <b>Clone Number</b>        | MCM6/3000                                                                                                                                                                        |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                          |                                                                                                                                         |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Expression System</b> | Bioreactor                                                                                                                              |
| <b>Storage</b>           | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>              | For research use only                                                                                                                   |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                                                                                         |

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)