

Product Datasheet

Mouse IgG2a FITC Negative Control Antibody (orb613769)

Catalog Number	orb613769
Category	Antibodies
Description	Mouse IgG2a FITC Negative Control Antibody
Species/Host	Mouse
Isotype	IgG2a
Conjugation	FITC
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide and 0.2% carrier protein
Concentration	See vial for concentration
Purification	Protein A/G Chromatography
Tested applications	FC
Dilution range	Titered for flow cytometry
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)