

## Product Datasheet

### BDNF Antibody (orb621597)

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb621597                                                                                                                               |
| <b>Category</b>            | Antibodies                                                                                                                              |
| <b>Description</b>         | BDNF Antibody                                                                                                                           |
| <b>Clonality</b>           | Recombinant                                                                                                                             |
| <b>Species/Host</b>        | Rabbit                                                                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Reactivity</b>          | Human, Mouse, Rat, Zebrafish                                                                                                            |
| <b>UniProt ID</b>          | <b>P23560</b>                                                                                                                           |
| <b>Tested applications</b> | ICC, IHC, WB                                                                                                                            |
| <b>Antibody Type</b>       | Primary Antibody                                                                                                                        |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)