

## Product Datasheet

### MOUSE ANTI SHEEP CD45R (orb623011)

|                             |                                                                                                                                         |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb623011                                                                                                                               |
| <b>Category</b>             | Antibodies                                                                                                                              |
| <b>Description</b>          | Mouse monoclonal antibody to CD45                                                                                                       |
| <b>Target</b>               | CD45R                                                                                                                                   |
| <b>Clonality</b>            | Monoclonal                                                                                                                              |
| <b>Species/Host</b>         | Mouse                                                                                                                                   |
| <b>Isotype</b>              | IgG1                                                                                                                                    |
| <b>Conjugation</b>          | Unconjugated                                                                                                                            |
| <b>Reactivity</b>           | Sheep                                                                                                                                   |
| <b>Form/Appearance</b>      | Purified IgG - liquid                                                                                                                   |
| <b>Concentration</b>        | IgG concentration 1.0 mg/ml                                                                                                             |
| <b>Buffer/Preservatives</b> | Purified IgG - liquid                                                                                                                   |
| <b>Purity</b>               | Purified                                                                                                                                |
| <b>Immunogen</b>            | Sheep thymocytes                                                                                                                        |
| <b>Tested applications</b>  | FC                                                                                                                                      |
| <b>Clone Number</b>         | 20.96                                                                                                                                   |
| <b>Storage</b>              | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                 | For research use only                                                                                                                   |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Expiration Date**

12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)