

Product Datasheet

CD68 Antibody (RPE) (orb623022)

Catalog Number	orb623022
Category	Antibodies
Description	Mouse monoclonal antibody against CD68 conjugated to RPE
Target	CD68
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG2b
Conjugation	RPE
Reactivity	Human
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Buffer/Preservatives	Preservative: 0.09% Sodium Azide (NaN ₃)1% Bovine Serum Albumin5% Sucrose; Buffer: Phosphate buffered saline
Immunogen	Phytohaemagglutinin stimulated peripheral blood mononuclear cells.
UniProt ID	P34810
Tested applications	FC
Specificity	CD68
Antibody Type	Primary Antibody
Clone Number	Y1/82A

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

This product is shipped at ambient temperature. Prior to reconstitution store at +4°C. After reconstitution store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light.

Note

For research use only

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)