

Product Datasheet

Recombinant Mouse IL-1 beta (orb623192)

Description	IL-1 beta (IL-1 β) is a member of the interleukin 1 family of cytokines. The IL-1 beta cytokine is produced by activated macrophages as a proprotein, which is proteolytically processed to its active form by caspase 1 (CASP1/ICE). This cytokine is an important mediator of the inflammatory response, and is involved in a variety of cellular activities, including cell proliferation, differentiation, and apoptosis. Mouse IL-1 beta Recombinant Protein is purified interleukin-1 beta cytokine produced in yeast.
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Application notes	The Mouse APRIL protein can be used in cell culture, such as an APRIL ELISA Standard, and as a Western Blot Control.
Protein Sequence	VPIRQLHYRLRDEQQKSLVLSDPYELKALHLNGQNINQQVIFSMSFVQGEPNSNDKIPVAL GLKGKNLYLSCVMKDGTPTLQLESVDPKQYPKKKMEKRFVFNKIEVKSKEFESAEFPNW YISTSQAEHKPVFLGNNSGQDIIDFTMESVSS (152)
Source	Yeast
MW	17.4 kDa
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)