

Product Datasheet

Recombinant Mouse VEGF-A (orb623194)

Catalog Number	orb623194
Category	Proteins
Description	VEGF-A is member of the family of vascular endothelial growth factor (VEGF) proteins, which stimulate vasculogenesis and angiogenesis to restore oxygen supply to tissues. VEGF-A is also a vasodilator and increases microvascular permeability and was originally referred to as vascular permeability factor (VPF). Mouse VEGF-A Recombinant Protein is purified vascular endothelial growth factor A produced in yeast.
Form/Appearance	Lyophilized
Protein Sequence	APTTEGEQKSHEVIKFMVDVYQRSYCRPIETLVDIFQEYPDEIEYIFKPSCVPLMRCAGCC NDEALECVPTSESNITMQIMRIKPHQSQHIGEMSFLQHSRCECRPKKDRTKPENHCEPCS ERRKHLFVQDPQTCKCCKNTDSRCKARQLELNERTCRCDKPRR (164)
MW	19.3 kDa
Application notes	The Mouse IL-2 protein can be used in cell culture, as an IL-2 ELISA Standard, and as a Western Blot Control.
Source	Yeast
Storage	-20°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558