

Product Datasheet

Recombinant Mouse IL-1 alpha (orb623229)

Catalog Number	orb623229
Category	Proteins
Description	IL-1 alpha (IL-1 α , IL-1F1) is a member of the interleukin 1 family of cytokines. IL-1 alpha is an inflammatory cytokine active in the initiation of the inflammatory reaction and in driving Th1 and Th17 inflammatory responses. Mouse IL-1 alpha Recombinant Protein is purified interleukin-1 alpha cytokine produced in yeast.
Form/Appearance	Lyophilized
Protein Sequence	SAPYTYQSDL RYKLMKLVQR KFMNDSLNLQ TIYQDVDKHY LSTTWLNDLQ QEVKFDMYAY SSGGDDSKYP VTLKISDSQL FVSAQGEDQP VLLKELPETP KLITGSETDL IFFWKSINSK NYFTSAAYPE LFIATKEQSR VHLARGLPSM TDFQIS (156)
MW	18.0 kDa
Application notes	The rat CXCL9/MIG protein can be used in cell culture, as a CXCL9/MIG ELISA Standard, and as a Western Blot Control.
Source	Yeast
Storage	-20°C
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)