

Product Datasheet

Anti-Ag85b (M. tuberculosis) Purified (orb623563)

Description	Rabbit polyclonal antibody to Ag85b (M. tuberculosis)
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mycobacterium tuberculosis Ag85b produced in E. coli
Target	Ag85b (M. tuberculosis)
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit polyclonal
Clonality	Monoclonal
Clone Number	PAb (189)
Antibody Type	Primary Antibody
Purity	Purified by ligand affinity chromatography.
Uniprot ID	I6YBV4
Entrez	885785
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)