

Product Datasheet

Flag tag Mouse Monoclonal Antibody (orb627074)

Catalog Number	orb627074
Category	Antibodies
Description	Mouse monoclonal antibody to Flag tag
Target	Flag tag
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1
Conjugation	Unconjugated
Form/Appearance	liquid
Purity	≥95% as determined by SDS-PAGE
Purification	Protein A+G purified
Immunogen	Flag tag
Tested applications	ELISA, IP, WB
Dilution range	WB: 1:500-1:2000; IP: 1:500-1:2000
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States


Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date 12 months from date of receipt.



IP Result of anti-Flag tag (IP: orb627074, 3 ug; Detection: orb627074 1:2000) with Recombinant protein protein lysate 500 ug.



Flagged fusion protein were subjected to SDS PAGE followed by western blot with orb627074 (Flag tag antibody) at dilution of 1:5000.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)