

## Product Datasheet

### SCAF8 Antibody (orb630562)


|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to RBM16                                                                                                     |
| <b>Species/Host</b>        | Rabbit                                                                                                                                  |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | ELISA, IF, IHC, WB                                                                                                                      |
| <b>Immunogen</b>           | RNA binding motif protein 16                                                                                                            |
| <b>Target</b>              | SCAF8                                                                                                                                   |
| <b>Form/Appearance</b>     | liquid                                                                                                                                  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Isotype</b>             | IgG                                                                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                                                                              |
| <b>Antibody Type</b>       | Primary Antibody                                                                                                                        |
| <b>Purity</b>              | ≥95% as determined by SDS-PAGE                                                                                                          |
| <b>MW</b>                  | 170 kDa                                                                                                                                 |
| <b>Uniprot ID</b>          | <b>Q9UPN6</b>                                                                                                                           |
| <b>Dilution Range</b>      | WB: 1:1000-1:10000; IP:1:200-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100                                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                         |

#### Biorbyt Ltd.


7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



HT-1080 cells were subjected to SDS PAGE followed by western blot with orb630562 (RBM16 antibody) at dilution of 1:2000.



Immunohistochemistry of paraffin-embedded human testis using orb630562 (RBM16 antibody) at dilution of 1:100.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)