

## Product Datasheet

### UXT Rabbit Polyclonal Antibody (orb632469)

|                            |                                   |
|----------------------------|-----------------------------------|
| <b>Catalog Number</b>      | orb632469                         |
| <b>Category</b>            | Antibodies                        |
| <b>Description</b>         | Rabbit polyclonal antibody to UXT |
| <b>Target</b>              | UXT                               |
| <b>Clonality</b>           | Polyclonal                        |
| <b>Species/Host</b>        | Rabbit                            |
| <b>Isotype</b>             | IgG                               |
| <b>Conjugation</b>         | Unconjugated                      |
| <b>Reactivity</b>          | Human, Mouse, Rat                 |
| <b>Form/Appearance</b>     | liquid                            |
| <b>Purity</b>              | ≥95% as determined by SDS-PAGE    |
| <b>Purification</b>        | Immunogen affinity purified       |
| <b>Immunogen</b>           | ubiquitously-expressed transcript |
| <b>UniProt ID</b>          | <b>Q9UBK9</b>                     |
| <b>MW</b>                  | 78 kDa, 36 kDa                    |
| <b>Tested applications</b> | ELISA, WB                         |
| <b>Dilution range</b>      | WB: 1:200-1:1000                  |
| <b>Antibody Type</b>       | Primary Antibody                  |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                        |   |
|------------------------|---|
| <b>Storage</b>         | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>            | For research use only   |
| <b>Expiration Date</b> | 12 months from date of receipt.   |



HeLa cells were subjected to SDS PAGE followed by western blot with orb632469 (UXT antibody) at dilution of 1:400.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)