

## Product Datasheet

### CD3E Rabbit Polyclonal Antibody (orb633084)

|                            |  |
|----------------------------|--|
| <b>Catalog Number</b>      | orb633084                                |
| <b>Category</b>            | Antibodies                               |
| <b>Description</b>         | Rabbit polyclonal antibody to CD3E       |
| <b>Target</b>              | CD3E                                     |
| <b>Clonality</b>           | Polyclonal                               |
| <b>Species/Host</b>        | Rabbit                                   |
| <b>Isotype</b>             | IgG                                      |
| <b>Conjugation</b>         | Unconjugated                             |
| <b>Reactivity</b>          | Human                                    |
| <b>Form/Appearance</b>     | liquid                                   |
| <b>Purity</b>              | ≥95% as determined by SDS-PAGE           |
| <b>Purification</b>        | Immunogen affinity purified              |
| <b>Immunogen</b>           | CD3e molecule, epsilon (CD3-TCR complex) |
| <b>UniProt ID</b>          | <b>P07766</b>                            |
| <b>MW</b>                  | 23 kDa                                   |
| <b>Tested applications</b> | ELISA, WB                                |
| <b>Dilution range</b>      | WB: 1:500 - 1:2000                       |
| <b>Antibody Type</b>       | Primary Antibody                         |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Storage**


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Expiration Date**

12 months from date of receipt.



various lysates were subjected to SDS PAGE followed by western blot with orb633084 (CD3E Antibody) at dilution of 1:1000.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)