

## Product Datasheet

### Streptavidin antibody (HRP) (orb633180)

<b>Catalog Number</b>	orb633180
<b>Description</b>	Rabbit polyclonal antibody against Streptavidin conjugated to HRP
<b>Reactivity</b>	Bacteria
<b>Tested Applications</b>	ELISA, WB
<b>Preservatives</b>	Lyophilized from a 0.2 µm filtered solution in PBS with 5% trehalose, pH 7.4
<b>Storage</b>	The lyophilized protein is stable at -20°C for up to 1 year, the protein solution can be stored at 2-8°C for up to 1 week under sterile condition. For extended storage, it is recommended to further dilute in a buffer containing 0.1% BSA (make sure BSA will NOT interfere with your experiment ) and store in working aliquots at -20°C to -80°C for 6 months. Avoid repeated freeze/thaw cycle.
<b>Tag</b>	N-terminal His Tag
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge the vial prior to opening, reconstitute in sterile distilled water or buffer of interest to concentration of 0.1-1mg/ml by gently pipetting 2-3 times, don't vortex
<b>Antibody Type</b>	Secondary Antibody
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE
<b>Source</b>	E.Coli
<b>Uniprot ID</b>	<b>P22629</b>
<b>Expiration Date</b>	6 months from date of receipt.

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)