

Product Datasheet

HRP conjugated Goat Rat IgG (H+L) Antibody (orb633194)

Catalog Number	orb633194
Category	Antibodies
Description	Goat polyclonal antibody against Goat Anti-Rat IgG (H+L) conjugated to HRP
Clonality	Polyclonal
Species/Host	Goat
Isotype	IgG
Conjugation	HRP
Reactivity	Rat
Form/Appearance	liquid
Concentration	3mg/ml BSA (IgG-Free, Protease-Free)
Buffer/Preservatives	3mg/ml rAlbumin (IgG-Free, Protease-Free)
Tested applications	ELISA, WB
Dilution range	1:2000-1:10000 for western blotting with ECL substrates1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates
Antibody Type	Secondary Antibody
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558