

Product Datasheet

Donkey Anti-Rat IgG (H+L) Secondary Antibody (AP) (orb633204)

Catalog Number	orb633204
Category	Antibodies
Description	Rat monoclonal antibody against Affinipure Donkey Anti- Rat IgG (H+L) conjugated to AP
Clonality	Polyclonal
Species/Host	Donkey
Isotype	IgG
Conjugation	AP
Reactivity	Rat
Form/Appearance	liquid
Buffer/Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Tested applications	ELISA, WB
Dilution range	1:1000-1:3000 for ELISA and Western blotting
Antibody Type	Secondary Antibody
Storage	0.01M Tris-HCl, 0.25M NaCl, 50% glycerol, pH8.0. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)