

Product Datasheet

Rabbit Anti-Goat IgG (H+L) Secondary Antibody (FITC) (orb633212)

Catalog Number	orb633212
Category	Antibodies
Description	Rabbit polyclonal antibody against Affinipure Rabbit Anti-Goat IgG (H+L) conjugated to FITC
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	FITC
Reactivity	Goat
Form/Appearance	liquid
Concentration	5mg/ml BSA (IgG-Free, Protease-Free)
Buffer/Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Tested applications	FC, IF
Dilution range	1:20-1:100 for most applications
Antibody Type	Secondary Antibody
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol□pH7.6. -20°C for 24 months.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)