

Product Datasheet

Biotin conjugated Affinipure Donkey Anti-Goat IgG (H+L) (orb633220)

Catalog Number	orb633220
Description	Goat polyclonal antibody against Affinipure Donkey Anti-Goat IgG (H+L) conjugated to Biotin
Species/Host	Donkey
Reactivity	Goat
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC, WB
Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Antibody Type	Secondary Antibody
Dilution Range	1:100 – 1:1000 for enzyme immunohistochemistry on tissue sections1:4000-1:80000 for ELISA and Western blotting with chromogenic substrates
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)