

Product Datasheet

Donkey Anti-Rabbit IgG (H+L) Secondary Antibody (Biotin) (orb633226)

Catalog Number	orb633226
Category	Antibodies
Description	Rabbit polyclonal antibody against Affinipure Donkey Anti- Rabbit IgG (H+L) conjugated to Biotin
Clonality	Polyclonal
Species/Host	Donkey
Isotype	IgG
Conjugation	Biotin
Reactivity	Rabbit
Form/Appearance	liquid
Concentration	5mg/ml BSA (IgG-Free, Protease-Free)
Buffer/Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Tested applications	ELISA, ICC, IHC, WB
Dilution range	1:100 - 1:1000 for enzyme immunohistochemistry on tissue sections 1:4000-1:80000 for ELISA and Western blotting with chromogenic substrates
Antibody Type	Secondary Antibody
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)