

Product Datasheet

AMCA conjugated Affinipure Donkey Goat IgG (H+L) Antibody (orb633256)

Catalog Number	orb633256
Category	Antibodies
Description	Goat polyclonal antibody against Affinipure Donkey Anti-Goat IgG (H+L) conjugated to AMCA
Clonality	Polyclonal
Species/Host	Donkey
Isotype	IgG
Conjugation	AMCA
Reactivity	Goat
Form/Appearance	liquid
Buffer/Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Tested applications	FC, IF
Dilution range	1:20-1:100
Antibody Type	Secondary Antibody
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)