

Product Datasheet

GRM2 Antibody (orb636229)

Description	Mouse monoclonal antibody to GRM2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Purified recombinant fragment of human GRM2 (AA: extra 414-558) expressed in E. Coli.
Target	GRM2
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Form/Appearance	Purified antibody in PBS with 0.05% sodium azide.
Concentration	p>1mg/ml
Storage	This product will ship in a box containing blue ice at a temperature of 4°C.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Clone Number	6E4F6
Antibody Type	Primary Antibody
MW	95.6kDa
Uniprot ID	Q14416

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Entrez


2912

Dilution Range


ELISA: 1/10000, FACS: 1/200 - 1/400, WB: 1/500 - 1/2000

Expiration Date


12 months from date of receipt.



Western blot analysis of SK-N-SH (1), Jurkat (2), and SW620 (3) cell lysate using GRM2 antibody.



Line graph illustrates about the Ag-Ab reactions using different concentrations of antigen and serial dilutions of GRM2 Antibody.




Flow cytometric analysis of SK-N-SH cells using GRM2 antibody.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)



Western blot analysis of recombinant protein using GRM2 antibody.



Western blot analysis of HEK293 (1) & HEK293 (2) cell lysate using GRM2 antibody.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)