

Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleocapsid Mouse Monoclonal Antibody (orb638871)

Catalog Number	orb638871
Category	Antibodies
Description	SARS-CoV-2 (2019-nCoV) Nucleocapsid Mouse Monoclonal Antibody
Target	S
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG2b
Conjugation	Unconjugated
Reactivity	Virus
Predicted Reactivity	Virus
Form/Appearance	Liquid
Concentration	≥1 mg/ml
Buffer/Preservatives	0.01M PBS (pH7.4) with 0.02% Proclin300.
Purification	Affinity purified by Protein G
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protein full length
MW	46 kDa
Tested applications	WB

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Dilution range	WB=1:500-2000
Antibody Type	Primary Antibody
Clone Number	7A7
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only



Sample: Lane 1: Recombinant SARS-CoV-2 N Protein (Omicron (B.1.1.529) (P13L, E31del, R32del, S33del, R203K, G204R) at 100 ng, Lane 2: Recombinant SARS-CoV-2 N Protein (Omicron (B.1.1.529) (P13L, E31del, R32del, S33del, R203K, G204R) at 50 ng, Lane 3: Recombinant SARS-CoV-2 N Protein (Omicron (B.1.1.529) (P13L, E31del, R32del, S33del, R203K, G204R) at 25 ng, Lane 4: Recombinant SARS-CoV-2 N Protein (Omicron (B.1.1.529) (P13L, E31del, R32del, S33del, R203K, G204R) at 10 ng, Lane 5: Recombinant SARS-CoV-2 N Protein (Omicron (B.1.1.529) (P13L, E31del, R32del, S33del, R203K, G204R) at 5 ng, Lane 6: Recombinant SARS-CoV-2 N Protein (WT) at 10 ng, Primary: Mouse Anti-SARS-CoV-2 N protein (orb638871) at 1/1000 dilution, Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution, Predicted band size: 45.3 kDa, Observed band size: 50 kDa.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Sample: Lane 1: Recombinant SARS-CoV-2 N protein (WT) (His Tag), Lane 2: Recombinant SARS-CoV-2 N protein (Q9H, P67S, P80R, P151L, S183Y) (His Tag), Lane 3: Recombinant SARS-CoV-2 N protein (D3L, P13T, D103Y, D128Y, H145Y, R203K, G204R, T205I, S235F) (His Tag), Lane 4: Recombinant SARS-CoV-2 N protein (Del204, Del215) (His Tag), Lane 5: Recombinant SARS-CoV-2 N protein (R203M, D377Y) (His Tag), Lane 6: Recombinant SARS-CoV-2 (Omicron, B.1.1.529) N protein (P13L, E31del, R32del, S33del, R203K, G204R) (N-His Tag), Primary: Anti-SARS-CoV-2 (2019-nCoV) Nucleocapsid (orb638871) at 1/1000 dilution, Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution, Predicted band size: 46 kDa, Observed band size: 50 kDa.



Sample: Lane 1: SARS-CoV-2 N protein (orb1516719) at 50 ng, Lane 2: SARS-CoV-2 N protein (orb1516719) at 25 ng, Lane 3: SARS-CoV-2 N protein (orb1516719) at 5 ng, Lane 4: SARS-CoV-2 N protein (orb1516719) at 2 ng, Lane 5: SARS-CoV-2 N protein (orb1516719) at 1 ng, Lane 6: SARS-CoV-2 N protein (orb1516719) at 0.5 ng, Primary: Anti-SARS-CoV-2 N protein at 1/1000 dilution, Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution, Predicted band size: 46 kD, Observed band size: 46 kD.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558