

Product Datasheet

NUP50 Antibody (ascites) (orb652075)

Catalog Number	orb652075
Category	Antibodies
Description	Mouse monoclonal antibody to NUP50
Target	NUP50
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1,k
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Mouse monoclonal antibody supplied in crude ascites with 0.09% (W/V) sodium azide.
Immunogen	Recombinant Protein
UniProt ID	Q9UKX7
MW	50144
Tested applications	WB
Dilution range	WB: 1:1000-3200
Specificity	This NUP50 monoclonal antibody is generated from mouse immunized with NUP50 recombinant protein.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Expiration Date	12 months from date of receipt.



NUP50 western blot analysis in Jurkat cell line lysates (35 µg/lane). This demonstrates the NUP50 antibody detected the NUP50 protein (arrow).

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)