


Product Datasheet

Cell Supernatant Exosome Isolation Kit (orb653766)

Catalog Number	orb653766
Category	Assays and Kits
Description	Exosomes are a type of extracellular vesicles secreted by cells that contain a variety of RNA, proteins, and phospholipids with a diameter of about 30-150 nm. Exosomes play an important role in the transfer of material and information between cells.
Sample Types	Cell Culture Supernatant, Serum, Plasma, Breast Milk, Urine, Saliva, Yeast, Plant
Storage	Storage: Store at 2-8°C. Shelf Life: One year from the date of receipt
Hazard Information	Non-Toxic
Note	For research use only
Expiration Date	12 months from date of receipt.



CD9 (orb688874); 1. Exosomes extracted from plasma; 2. Exosomes extracted from serum; 3. Exosomes extracted from urine; 4. Exosomes extracted from Hela cells; 5. Exosomes extracted from Breast Milk; 6. Exosomes extracted from saliva;

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)