

Product Datasheet

JAG1 (N) Antibody, Rabbit Polyclonal (orb66697)

Catalog Number	orb66697
Category	Antibodies
Description	Rabbit polyclonal to JAG1. Jagged 1 (JAG1) protein is the ligand for the receptor Notch 1. JAG1 signaling through Notch 1. Notch signaling is important for embryonic development, hematopoiesis, and stem cell maintenance. Mutations that alter JAG1 have been implicated in Alagille syndrome type 1.
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human, Mouse
Buffer/Preservatives	0.01% NaN ₃
Purity	Affinity purification
Immunogen	N-terminal region of human JAG1
UniProt ID	P78504
MW	134
Tested applications	WB
Dilution range	WB: 1:1000-1:5000
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Western blot analysis of total cell extracts from mouse liver using JAG1 antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)