

## Product Datasheet

### Dexas1 Rabbit Polyclonal Antibody (orb667895)

|                             |   |
|-----------------------------|---|
| <b>Catalog Number</b>       | orb667895   |
| <b>Category</b>             | Antibodies  |
| <b>Description</b>          | Rabbit polyclonal antibody to Dexas1  |
| <b>Target</b>               | RASD1   |
| <b>Clonality</b>            | Polyclonal  |
| <b>Species/Host</b>         | Rabbit  |
| <b>Conjugation</b>          | Unconjugated  |
| <b>Reactivity</b>           | Human, Mouse  |
| <b>Form/Appearance</b>      | Liquid  |
| <b>Buffer/Preservatives</b> | 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| <b>Purification</b>         | The antibody was purified by immunogen affinity chromatography.                                 |
| <b>Immunogen</b>            | Recombinant fusion protein of human Dexas1. The exact sequence is proprietary.                  |
| <b>UniProt ID</b>           | <b>Q9Y272, O35626</b>   |
| <b>Tested applications</b>  | WB  |
| <b>Dilution range</b>       | WB: 1-500-2000  |
| <b>Specificity</b>          | Recognizes endogenous levels of Dexas1 protein  |
| <b>Antibody Type</b>        | Primary Antibody  |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                        |   |
|------------------------|---|
| <b>Source</b>          | Rabbit  |
| <b>Storage</b>         | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>            | For research use only   |
| <b>Entrez</b>          | <b>51655, 19416</b>   |
| <b>Expiration Date</b> | 12 months from date of receipt.   |



Western blot analysis of Dexras1 expression in HEK293T (A), mouse brain (B), mouse heart (C) whole cell lysates. (Predicted band size: 13; 31 kD; Observed band size: 31 kD)

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)